

Supplementary information on EN 18032 (informative)

Foodstuff — Quick Method for the Analysis of Multiple Highly Polar Pesticides and their Metabolites in Foodstuff Involving Extraction with Acidified Methanol and Measurement by LC- or IC-MS/MS (QuPPE-Method) - Supplementary information on the method

Lebensmittel — Schnellmethode zur Bestimmung mehrerer hochpolarer Pestizide und ihrer Metaboliten in Lebensmitteln nach Extraktion mit angesäuertem Methanol und Messung mittels LC- oder IC-MS/MS (QuPPE-Methode) - Ergänzende Informationen zur Methode


Produit alimentaire — Méthode rapide pour l'analyse de plusieurs pesticides hautement polaires et de leurs métabolites dans les aliments impliquant une extraction avec du méthanol acidifié et une mesure par LC- ou IC-MS/MS (QuPPE-Method) - Informations complémentaires sur la méthode

Contents	Page
Supplement A (informative): Schemes showing the method at a glance	3
SA.1 QuPpe Method at a glance (for most fruits and vegetables)	4
SA.2 QuPpe Method at a glance (for E4 honey)	5
SA.3 QuPpe Method at a glance (for E5b cereals, pulses, nuts and oily seeds)	6
SA.4 QuPpe Method at a glance (for E6a liver, milk, kidney, egg and muscle)	7
SA.5 QuPpe Method at a glance (for E6b liver, milk, kidney, egg and muscle)	8
SA.6 QuPpe Method at a glance (for E7 animal fat)	9
SA.7 Decision trees on the use of quantification modules	10
Supplement B (informative): Validation Data	12
Supplement C (informative): Additional information	13
Table SC.1 — Conversion factors between typical purchased standards and target analytes	14
Table SC.2 — Exemplary concentrations and solvents of stock and working solutions of analytes and their respective IL-ISs	15


Supplement A (informative):

Schemes showing the method at a glance


SA.1 QuPPE Method at a glance (for most fruits and vegetables)




SA.2 QuPPE Method at a glance (for E4 honey)




SA.3 QuPPE Method at a glance (for E5b cereals, pulses, nuts and oily seeds)




SA.4 QuPPE Method at a glance (for E6a liver, milk, kidney, egg and muscle)



SA.5 QuPPE Method at a glance (for E6b liver, milk, kidney, egg and muscle)




SA.6 QuPPE Method at a glance (for E7 animal fat)



SA.7 Decision trees on the use of quantification modules


Decision tree
where (relative) matrix effects are NOT negligible including after dilution



- Q1a: Solvent-based without IL-IS
 - Q1b: Solvent-based with IL-IS added to the test portion
 - Q2a: Matrix-based without IL-IS added to the test portion
 - Q2b: Matrix-based with IL-IS added to the test portion
 - Q3a: Matrix-matched without IL-IS
 - Q3b: Matrix-matched with IL-IS added to the test portion
 - Q4a: Procedural calibration without IL-IS added to the test portion or without IL-IS but in a matching matrix
 - Q4b: Procedural calibration with IL-IS added to the test portion or without IL-IS but in a matching matrix
 - Q5a: Standard addition on extract aliquots, optionally without IL-IS added to the test portion
 - Q5b: Standard addition on extract aliquots, optionally with IL-IS added to the test portion
 - Q6a: Standard addition on sample portions, optionally without IL-IS added to the test portion
 - Q6b: Standard addition on sample portions, optionally with IL-IS added to the test portion
- *In all these cases IL-IS may also be added to an aliquot of the extracts (but water adjustment is necessary)
 ** 2nd choice approach is in some cases preferable to the first one (e.g. less work-intensive or easier to handle in routine, e.g. matrix-based calibrations more practical than matrix-matched ones, within a sequence of mixed commodities)

		Options:				
		preferred	second**	also possible		
			Q3b	Q4b	Q5b	Q5b
			Q3a	Q4a	Q5a	Q6a
Q1b		Q2b		Q4b	Q5b	Q6b
				Q5a	Q6a	
			Q3b	Q4b		Q6b
				Q4a		Q6a
Q1b		Q2b		Q4b		Q6b
						Q6a

miro



Q1a: Solvent-based without IL-IS
Q1b: Solvent-based with IL-IS added to the test portion
Q2a: Matrix-based without IL-IS added to the test portion
Q2b: Matrix-based with IL-IS added to the test portion
Q3a: Matrix-matched without IL-IS
Q3b: Matrix-matched with IL-IS added to the test portion
Q4a: Procedural calibration without IL-IS added to the test portion or without IL-IS but in a matching matrix
Q4b: Procedural calibration with IL-IS added to the test portion or without IL-IS but in a matching matrix
Q5a: Standard addition on extract aliquots, optionally without IL-IS added to the test portion
Q5b: Standard addition on extract aliquots, optionally with IL-IS added to the test portion
Q6a: Standard addition on sample portions, optionally without IL-IS added to the test portion
Q6b: Standard addition on sample portions, optionally with IL-IS added to the test portion
 *In all these cases IL-IS may also be added to an aliquot of the extracts (but water adjustment is necessary)
 **2nd choice approach is in some cases preferable to the first one (e.g. less work-intensive or easier to handle in routine, e.g. matrix-based calibrations more practical than matrix-matched ones, within a sequence of mixed commodities)

		Options:					
		preferred	second**	also possible			
				Q3b	Q4b	Q5b	Q6b
			Q3a	Q4a	Q5a	Q6a	
	Q1b	Q2b		Q4b	Q5b	Q6b	
Q1a	Q2a			Q4a	Q5a	Q6a	
			Q3b	Q4b			Q6b
				Q4a		Q6a	
	Q1b	Q2b		Q4b			Q6b
				Q4a		Q6a	

miro

Supplement B (informative):

Validation Data

Validation data from interlaboratory validation studies round 1 to 3 is available in tabular form as
Excel-Sheet here: [QuPPE.eu - EN Supplements](https://quppe.eu/en/supplements)

Supplement C (informative):

Additional information

Table SC.1 — Conversion factors between typical purchased standards and target analytes

Compound	MW [g/mol]	Compound as sold	MW [g/mol]	Conversion Factor (CF)	Inverse CF
Bromide (anion)	79,9	Potassium bromide	119,0	0,67	1,49
Chlorate (anion)	83,5	Chlorate-sodium	106,4	0,78	1,27
Chlormequat (cation) ^a	122,6	Chlormequat-chloride ^a	158,1	0,78	1,29
Chlormequat-D ₄ (cation)	126,6	Chlormequat-D ₄ -chloride	162,1	0,78	1,28
Difenzoquat (cation)	249,3	Difenzoquat-methylsulfate	360,4	0,69	1,45
Difluoroacetic acid- ¹³ C ₂	96,0	Sodium difluoroacetate- ¹³ C ₂	120,0	0,80	1,25
Fosetyl	110,0	Fosetyl-Al	118,0 ^b	0,93	1,07
Fosetyl-D ₅	115,0	Fosetyl-D ₅ -Al	123,0 ^b	0,93	1,07
		Fosetyl-D ₅ -sodium	137,0	0,84	1,19
Glufosinate	181,1	Glufosinate-ammonium	198,2	0,91	1,09
Glufosinate-D ₃	184,1	Glufosinate-D ₃ -ammonium hydrate	243,2	0,76	1,32
		Glufosinate-D ₃ -hydrochloride	220,6	0,83	1,20
Mepiquat (cation) ^a	114,2	Mepiquat-chloride ^a	149,7	0,76	1,31
Mepiquat-D ₃ (cation)	117,2	Mepiquat-D ₃ -iodide	244,1	0,48	2,08
Mepiquat-4-hydroxy	130,2	Mepiquat-4-hydroxy-chloride	165,7	0,79	1,27
N-Acetyl-Glufosinate	223,2	N-Acetyl-Glufosinate-disodium	267,2	0,84	1,20
N-Acetyl-Glufosinate-D ₃	226,2	N-Acetyl-Glufosinate-D ₃ -disodium	270,2	0,84	1,19
Nereistoxin	149,3	Nereistoxin-oxalate	239,3	0,62	1,60
Nereistoxin-D ₆	155,3	Nereistoxin-D ₆ -oxalate	245,3	0,63	1,58
Nicotine	162,2	Nicotine hemisulfate	422,5 ^c	0,77	1,30
Propamocarb-N-oxide	204,3	Propamocarb-N-oxide hydrochloride	240,7	0,85	1,17
Trifluoroacetic acid	113,0	Sodium trifluoroacetate	136,0	0,83	1,20
Trifluoroacetic acid- ¹³ C ₂	115,0	Sodium trifluoroacetate- ¹³ C ₂	138,0	0,83	1,20
Trimethylsulfonium (cation)	77,2	Trimethylsulfonium-iodide	204,1	0,38	2,64
Trimethylsulfonium-D ₉ (cation)	86,2	Trimethylsulfonium-D ₉ -iodide	213,1	0,40	2,47

^a Attention: The current EU – MRLs are expressed as the respective chloride salts. Thus, no mathematical conversion of the chloride to the cation is needed.

^b The MW given here includes only 1/3 of the aluminium mass (one aluminium ion is bound to three fosetyl anions)

^c The MW refers to the following formula (C₁₀H₁₄N₂)₂ · H₂SO₄ which entails two nicotine molecules.

Table SC.2 — Exemplary concentrations and solvents of stock and working solutions of analytes and their respective IL-ISs

Compound	Stock Solution (exemplary)		Working Solutions including mixtures (exemplary)	
	Solvent used to prepare	[mg/mL]	Solvent used to prepare	[µg/mL]
Aminocyclopyrachlor	MeOH	1	MeOH	10 / 1 / 0,1
Amitrole	MeOH	1	MeOH	10 / 1 / 0,1
AMPA	10 % ACN in water	1	10 % ACN in water	10 / 1 / 0,1
Bromide	MeOH	1	MeOH	10 / 1 / 0,1 / 0,01
Chlorate	MeOH	1	MeOH	10 / 1 / 0,1 / 0,01
Chloridazon-desphenyl	MeOH	1	MeOH	10 / 1 / 0,1
Chlormequat	MeOH	1	MeOH	10 / 1 / 0,1
Cyanuric acid	MeOH	1	10 % ACN in water	10 / 1 / 0,1
Cyromazine	MeOH	1	MeOH	10 / 1 / 0,1
Daminozide	MeOH	1	MeOH	10 / 1 / 0,1
Difenzoquat	ACN	1	MeOH	10 / 1 / 0,1
Difluoroacetic acid	ACN with 5% water	1	ACN with 5% water	10 / 1 / 0,1
Ethephon	10 % ACN in water + 0,1 % HCl	1	10 % ACN in water + 0,1 % HCl	10 / 1 / 0,1
ETU	MeOH	1	MeOH	10 / 1 / 0,1
Fosetyl	10 % ACN in water	0,1	10 % ACN in water	10 / 1 / 0,1
Glufosinate	10 % ACN in water	1	10 % ACN in water	10 / 1 / 0,1
Glyphosate ^b	10 % ACN in water	1	10 % ACN in water	10 / 1 / 0,1
HEPA	10 % ACN in water	1	10 % ACN in water	10 / 1 / 0,1
Matrine	ACN	1	ACN	10 / 1 / 0,1
Maleic Hydrazide	MeOH	1	10 % ACN in water	10 / 1 / 0,1
Melamine	MeOH:water (90:10)	1	MeOH	10 / 1 / 0,1
Mepiquat	MeOH	1	MeOH	10 / 1 / 0,1
Mepiquat-4-hydroxy	MeOH	1	MeOH	10 / 1 / 0,1
MPPA	10 % ACN in water	1	10 % ACN in water	10 / 1 / 0,1
N-Acetyl-Glufosinate	10 % ACN in water	1	10 % ACN in water	10 / 1 / 0,1
N-Acetyl-Glyphosate	10 % ACN in water	1	10 % ACN in water	10 / 1 / 0,1
Nereistoxin	MeOH / water (3:1)	1	MeOH	10 / 1 / 0,1
Nicotine ^b	ACN	1	ACN	1 / 0,1
Oxymatrine	ACN	1	ACN	10 / 1 / 0,1

Supplementary information on EN 18032

Compound	Stock Solution (exemplary)		Working Solutions including mixtures (exemplary)	
	Solvent used to prepare	[mg/mL]	Solvent used to prepare	[µg/mL]
Perchlorate	MeOH	1	MeOH	10 / 1 / 0,1 / 0,01
Phosphonic acid ^b	Water (¹⁸ O-H ₂ O for the ILIS)	1	ACN ^c	10 / 1 / 0,1 / 0,01
Propamocarb	ACN	1	MeOH	10 / 1 / 0,1
Propamocarb-N-desmethyl	ACN:Acetone (1 mL acetone to initially dissolve)	1	MeOH	10 / 1 / 0,1
Propamocarb-N-oxide	MeOH	1	MeOH	10 / 1 / 0,1
PTU	MeOH	1	MeOH	10 / 1 / 0,1
Triazole-lactic acid	MeOH	1	MeOH	10 / 1 / 0,1
Triazole-acetic acid	MeOH	1	MeOH	10 / 1 / 0,1
Triazole-alanine	MeOH / water (1:3)	1	MeOH	10 / 1 / 0,1
Trifluoroacetic acid	ACN with 5% water	1	ACN with 5% water	10 / 1 / 0,1
Trimethylsulfonium	MeOH	1	MeOH	10 / 1 / 0,1
<p>^a Use plastic vials and protect solutions from light exposure.</p> <p>^b Use plastic vessels and stoppers for compounds that tend to interact with glass surfaces.</p> <p>^c Pure water (¹⁸O-H₂O for the IL-IS) is also suitable for the working solution. 10% ACN will reduce growth of microorganisms.</p> <p>MeOH: Methanol; ACN: Acetonitrile; FA: Formic acid.</p>				